Technovit® LC Xpress

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GB Instructions for use Technovit LC Xpress

Product description

Light-curing resin for creating sample blocks for material testing. For embedding all types of material samples that subsequently have to be ground and/or polished. Technovit LC Xpress works in combination with the Technotray POWER or Technotray POWER LED light-curing unit. The product is a chemical formulation. For the safe use of Technovit LC Xpress, the instructions on the packaging and/or the safety datasheet must be observed. Application in the laboratory by trained specialist personnel.

Properties

Technovit LC Xpress liquid is light-curing and is cured with non-hazardous, visible blue light in the Kulzer Technotray POWER or Technotray POWER LED polymerisation unit. The 1-component formulation means there is no mixing process and therefore no loss of material. As Technovit LC Xpress does not outgas, there is no need for a fume cupboard. Technovit LC Xpress produces orange-transparent, bubble-free specimens and can be optimally processed mechanically.

Preparation for embedding

Connect the light unit to the power supply and switch it on. Prepare the embedding moulds. The specimens to be embedded with Technovit LC Xpress must be clean and free of grease to prevent reactions with the embedding resin and to ensure optimum embedding. The long pot life of Technovit LC Xpress enables precise positioning of the sample, as polymerisation only starts when the material is irradiated with the corresponding light sources. Position the sample on the base of a Kulzer embedding mould corresponding to the sample size and, if necessary, secure it with Technovit LC Fix. Place the cylinder of the embedding mould on the base. Fill the mould up to approximately 3 mm below the edge with Technovit LC Xpress. The specimen must be completely covered with the embedding resin. To avoid air pockets under the sample, pour a level of Technovit LC Xpress into the PE or silicone mould to a height of approx. 3 mm. Place the sample in the liquid level and move it slightly back and forth. Place the sample in the centre and completely fill the mould.

Embedding with infiltration

In the case of porous samples or difficult geometric shapes of the sample, it may be advantageous to evacuate the prepared samples in a vacuum impregnation device without bubbles so that possible air inclusions can escape completely. For vacuum-assisted infiltration, only pour Technovit LC Xpress onto the relevant sample area. Then place the embedding mould together with the sample in the vacuum set and infiltrate. After infiltration is complete, fill the embedding mould with Technovit LC Xpress up to approximately 3 mm below the edge.

Polymerization

Insert the embedding moulds filled with Technovit LC Xpress into the corresponding light unit and start the device.

Working with Technotray POWER

Replace the two centre blue light tubes (top and bottom) with LED strips (item no. 66093673 “LED strip for Technotray POWER”). Switch on the device and open the drawer. Place the prepared moulds on the glass plate of the drawer. Set the time to “5 min” and close the drawer, the polymerisation process starts. At the end of the time programme, open the drawer and remove the finished samples.

Working with POWER LED

This device, which is compatible with Technovit LC Xpress, should be used wherever possible. Switch on the device and open the drawer. Place the prepared moulds on the glass plate (the device holds enough samples to fill a complete sample holder). Close the drawer. Select the “Technovit LC Xpress” programme setting on the device and start polymerisation. After the 90 second programme has elapsed, remove the samples.

Note: Depending on the volume of the material, samples may still be warm immediately after polymerisation – take care when removing them from the device!

Further processing

Remove the specimen blocks from the embedding mould and allow to cool further if necessary. Process further in accordance with the preparation protocol for the corresponding materials.

Delivery units

Technovit LC Xpress	– 1000 ml Liquid, unbreakable bottle
Technovit LC Fix	– 1 x 4 g Fixing paste
Light curing devices Technotray Power/Technotray POWER LED	
Transparent KULZER Embedding mould (PE or Silicone)	

Danger/Safety advices

Please follow the instructions on the product packaging and safety data sheets.

www.kulzer-technik.com

Our technical application advice, whether verbal, in writing or by means of tests, is given to the best of our knowledge, but is only non-binding information, including with regard to any third-party property rights, and does not release you from your obligation to test the products supplied by us for their suitability for the intended processes and purposes. We have no means of controlling the application, use or processing of the products, which are therefore the sole responsibility of the Buyer. It goes without saying that we guarantee the perfect quality of our products in accordance with our General Terms and Conditions of Sale and Delivery. If there is reason to assume that there are quality defects caused by the manufacturer, please inform us immediately. See the product packaging for the shelf life. Please comply with the temperatures specified in the processing instructions to ensure that the product is used correctly. The material must not be used after the expiry date. Do not store at a room temperature above 25 °C (77 °F). Protect from sunlight. Store in a dry place. Dispose of in accordance with official regulations. Do not dispose of contents and packaging that is not completely empty with household waste or into the sewage system. European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances. Surfaces or objects contaminated with Technovit LC Xpress can be cleaned with ethanol. Transport closed in the original container.

* For more information please visit www.kulzer-technik.com

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® = marca registrada, Kulzer GmbH, Hanau
Kulzer GmbH